



Technical specification

Width	220 mm
Length	210 mm
Thickness	51 mm
Weight	1450 g (battery pack included)

Sensors



miniflowmeter (code 900595)
for reusable and disposable turbine
dimension (Ø 30 mm, 42 mm)



Reusable soft, adult, MIR sensor for oximetry tests (code 919024) only for spirolab code 911081

Power supply	Rechargeable battery and mains power Ni-MH, 6 elements
Current capacity	4500 mAh
Consumption	average 250 mA
Backup battery voltage	none
Batteries charger	Output voltage=12 V, current=1A, compliant with EN 60601-1

Autonomy	~10 hours
Connectivity	USB 2.0, Bluetooth® 2.1
Display	7 inch colour touch screen LCD Display with 800x480 resolution

Keyboard	absent, touchscreen
Mouthpieces	Ø 30 mm (1.18 inch)
Type of electrical protection	Class II device
Safety level for shock hazard	Type BF Apparatus
Conditions of use	Apparatus for continuous use
Storage conditions	Temperature: MIN -40 °C, MAX +70 °C Humidity: MIN 10% RH; MAX 95%RH

Transport conditions	Temperature: MIN -40 °C, MAX +70 °C Humidity: MIN 10% RH; MAX 95%RH
-----------------------------	--

Operating conditions	Temperature: MIN + 10 °C, MAX + 40 °C Humidity: MIN 10% RH, MAX 95%RH
-----------------------------	--

Applied norms	Electrical Safety EN 60601-1 Electro Magnetic Compatibility EN 60601-1-2
----------------------	---

Degree of protection against water penetration	IPX1 appliance protected against water leaks
---	--

Codes and equipments

911080E0	spiro
911080E1	spiro with reusable turbine

911080E2	spiro with 120 FlowMir
911081E0	spiro+oxy
911081E1	spiro+oxy with reusable turbine
911081E2	spiro+oxy with 120 FlowMir

Spirometry

Flow sensor	bi-directional digital turbine
Volume rate	10 L
Flow range	±16L/s
Volume accuracy	±3% or 50 mL
Flow accuracy	±5% or 200 mL/s
Dynamic resistance	<0.5 cm H ₂ O/L/s
Temperature sensor	semiconductor (0-45°C)
Test available	FVC, VC, IVC, MVV, PRE-POST
Measured parameters	FVC, FEV1, FEV1/FVC%, FEV1/PEF, FEV1/VC, FEV1/FEF0.5, DTPEF, FEV 0.5, FEV0.5/FVC, FEV0.75, FEV0.75/FVC, FEV2, FEV2/FVC, FEV3, FEV3/FVC, FEV6, FEV1/FEV6, PEF, FEF25, FEF50, FEF75, FEF2575, FEF7585, FET, Vext, ELA, EVOL, FIVC, FIV1, PIF, FIV1/FIVC, FIF25, FIF50, FIF75, R50, MVVcal, PIF, IRV, VC, EVC, IVC, IC, ERV, IRV, FEV1/VC, TV, VE, RR, ti, te, ti/t-tot, TV/ti, MVV

Memory capacity

Up to 10000 tests

Oximetry (on request)

Measurement method	Red and infrared absorption
SpO2 range	0-99%
SpO2 accuracy	± 2% between 70-99% SpO2
Average number of heart beats for the %SpO2 calculation	8 beats
Pulse Rate range	18-300 BPM
Pulse Rate accuracy	± 2BPM or 2% whichever is greater
Average interval for the calculation of cardiac pulse	8 seconds

Signal quality indication	0 - 8 segments on display
Test available	spot
Measured parameters	SpO2% min, max, average BPM min, max, average Test duration % Bradycardia Duration (<40 BPM) % Tachycardia Duration (>120 BPM) % of Time with SpO2 ≤ 90% (T90%, T89%), T5

Memory capacity	about 500 hours oximetry
------------------------	--------------------------

Certificates & Registrations

CE 0476	MED 9826
FDA 510 (k)	K 052140
Health Canada	71191 (class II)
CND code	Z12150102 (spiro) Z1203020408 (spiro + oxy)
GMDN code	46906 (spiro), 45607 (spiro + oxy)
Ministry of Health	1272475/R (spiro) 1272476/R (spiro + oxy) 1645455/R (spiro)